



REF: 0807612

ANTIDECUBITUS CUSHION WITH GEL
ΜΑΞΙΛΑΡΙ ΚΑΤΑΚΛΙΣΕΩΝ GEL



www.mobiakcare.com



REF: 0807612

ANTIDECUBITUS CUSHION WITH GEL ΜΑΞΙΛΑΡΙ ΚΑΤΑΚΛΙΣΕΩΝ GEL

45 x 40 x 7cm



EN

- ✓ Memory Foam with Integrated Gel
- ✓ Cover with Zip

GR

- ✓ Memory Foam με Ενσωματωμένο Gel
- ✓ Κάλυμμα με Φερμουάρ

COMPOSITION / ΣΥΝΘΕΣΗ

INSIDE / ΕΣΩΤΕΡΙΚΟ : Regular Foam & Gel Pad

COVER / ΘΗΚΗ: PES 100%



www.mobiakcare.com

EN

The Washing Instructions Refer to the Cover ONLY:

- Hand Wash Only, using mild soap
- Do Not Bleach
- Do Not Iron
- Do Not Tumble Dry
- Do Not Dry Clean
- Do Not Dry Wring


CAUTION! THE CUSHION IS NOT WASHABLE

GR

Οδηγίες Πλυσίματος Αφορούν μόνο το Κάλυμμα:

- Πλύσιμο μόνο στο χέρι, χρησιμοποιώντας ήπιο απορρυπαντικό
- Απαγορεύεται το χλώριο
- Απαγορεύεται το σιδέρωμα
- Μην χρησιμοποιείτε στεγνωτήριο
- Απαγορεύεται το στεγνό καθάρισμα
- Απαγορεύεται το Στύψιμο

ΠΡΟΣΟΧΗ! ΤΟ ΜΑΞΙΛΑΡΙ ΔΕΝ ΠΛΕΝΕΤΑΙ



REF: 0807612

ANTIDECUBITUS CUSHION WITH GEL
ΜΑΞΙΛΑΡΙ ΚΑΤΑΚΛΙΣΕΩΝ GEL



LOT

MOBIAK S.A.
KATHIANA AKROTIRIOU CHANIA CRETE
GREECE / T: +30 28210 63222
e-mail: mobiakcaresales@mobiak.com
website: www.mobiakcare.com



REF: 0807612

ANTIDECUBITUS CUSHION WITH GEL
ΜΑΞΙΛΑΡΙ ΚΑΤΑΚΛΙΣΕΩΝ GEL



LOT

MOBIAK S.A.
KATHIANA AKROTIRIOU CHANIA CRETE
GREECE / T: +30 28210 63222
e-mail: mobiakcaresales@mobiak.com
website: www.mobiakcare.com

