



Heat and Moisture Exchange Filter

Instruction for Use 【Structure Performance】



STERILE EO



CE 0123



This product is mainly composed of filtering medium, such as medical polymeric shell, filtration membrane, etc.. After Ethylene Oxide Sterilization, the product is aseptic, non-toxic, and can be used safely and effectively.

【Storage Conditions】

Stored in a place whose relative humidity is not greater than 65%, with dry air, and no corrosive gas.

REF

ZF-001

LOT

25021501



20250215



20280214

【Points for Attention】

1. If the inner packing is broken, stop using;
2. The product is a one-time use product, repeated use is strictly prohibited;

Company: SHANDONG ZHENFU MEDICAL
DEVICE CO., LTD

Address: High-Tech Industrial Park, 276800
Rizhao City, P.R.C

TEL: +86-633-8268868 Fax: +86-633-8268988
Email: sales@zhenfumedical.com

EC REF

Company: Shanghai International
Holding Corp. GmbH (Europe)

Address: Eiffelstrasse 80, 20537
Hamburg, Germany

TEL: +49-40-2513175 Fax: +49-40-255726
Email: shholding@hotmail.com